


PROMO Product Information

Corrects IM angle and rotational deformity through a single, guided, oblique osteotomy


The PROMO system was developed based on the understanding that up to 87% of hallux valgus deformities have a frontal plane deformity¹

FEATURES & BENEFITS


- **Simultaneous correction of IM angle and rotational deformity**
 - Rotation in hallux valgus has been identified as a risk factor for recurrence²
- **Single oblique osteotomy with no wedge resection**
 - Joint sparing procedure with powerful correction
- **Guided jig system helps construct the osteotomy**
 - Repeatable and streamlined technique



Transverse Plane Correction



Frontal Plane Correction




PROMO Plate and Crossing Screw Fixation

THE PROMO SOLUTION


Positioning Jig

Positions K-wires to establish obliquity of the osteotomy




Cutting Jig

Guides osteotomy for patient specific correction



Rotation Jig

Indicates optimal rotational correction



Angle Values Table

- **Measurements needed:**
 - IM Angle
 - Rotation Angle
- **Calculations for the osteotomy cut angle are simplified into angle values table**
- **Values on table associated with jigs to guide construction of the osteotomy**

		Rotation Angle (°)			Osteotomy Cut Angle
		10-19	20-29	30-39	
IM Angle (°)	8-10	38	28	23	}
	11-12	47	33	28	
	13-14	55	38	33	
	15-17	55	42	38	
	18-20	55	47	42	

PROMO SPECIFIC PLATES

Straight



23-28-33-38



42-47-55

Angled



23-28-33-38



42-47-55

* Plate numbers correspond to osteotomy cut angles

Paragon 28, Inc.
4B Inverness Ct. E., Suite 280
Englewood, CO 80112 USA
(855) 786-2828

Paragon 28 Medical Devices Trading Limited
43 Fitzwilliam Square West
Dublin 2, D02 K792, Ireland
+353(0) 15414756


PROMO-01 RevC[™] Trademarks and [®] Registered Trademarks of Paragon 28, Inc. © Copyright 2018 Paragon 28, Inc. All rights reserved. Patents: www.paragon28.com/patents

1. Kim Y, Kim JS, Young KW, et al. A new measure of tibial sesamoid position in hallux valgus in relation to the coronal rotation of the first metatarsal in CT scans. Foot Ankle Int. 2015; 36(8): 944-52

2. Okuda R, Kinoshita M, et al. The shape of the lateral edge of the first metatarsal head as a risk factor for recurrence of hallux valgus. JBJS. 2007; 89(10):2163-72

PROMO Fixation

- **Medial wall plate**
 - Plate position allows for improved resistance to bending forces
- **Patent Pending Precision[®] Guide PROMO**
 - Allows for ideal cross screw placement while avoiding hardware collision
- **PROMO plate options designed around the length and angle of the osteotomy**



Precision[®] Guide
PROMO